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
1. HPMC (Hydroxy Propyl Methyl Cellulose)

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Concentration - viscosity relationship of PHARMACOAT® in water/ethanol (50:50)

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Relationship of the mixing ratio in the methylene chloride/ethanol system and the viscosity of 10% solutions of PHARMACOAT®

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Concentration - viscosity relationship of PHARMACOAT® in methylene chloride/ethanol (50:50)

1. HPMC (Hydroxy Propyl Methyl Cellulose)

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The effect of TiO₂ on the tensile strength of PHARMACOAT® films

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The effect of TiO₂ on the elongation of PHARMACOAT® films

1. HPMC (Hydroxy Propyl Methyl Cellulose)

Outline of types of Hypromellose (hydroxypropyl methylcellulose) or Methylcellulose

Labeled viscosity* ¹ (mPa·s)	Hypromellose Substitution type: 2208	Hypromellose Substitution type: 2906	Hypromellose Substitution type: 2910	Methylcellulose
3				
4				MCE - 4
4.5				
6			SE - 06	
15			SE - 15	MCE - 15
25				MCE - 25
50		SFE - 50	SE - 50	
100	NE - 100			MCE - 100
400		SFE - 400		MCE - 400
1500				MCE - 1500
4000	NE - 4000	SFE - 4000	SE - 4000	MCE - 4000
10000				
15000				
100000				
Substituent content Methoxy group Hydroxypropoxy group	19.0 - 24.0 %* ² 4.0 - 12.0 %* ²	27.0 - 30.0 % 4.0 - 7.5 %	28.0 - 30.0 % 7.0 - 12.0 %	26.0 - 33.0 % -

*¹ The viscosity is measured with 2 wt % aqueous solution at 20 °C.

Application Guide for HPMC

Application	Purpose
Film coating for tablet	Taste masking, Color masking, Coloring, Hardness (improve friability), Stability
Sub-coating for sugar coating	Stability, Prevention of moisture penetration to the core
Sub-coating for enteric coating	Prevention of interaction
Binder for tablet and granule	Compressibility, Compatibility

1. HPMC (Hydroxy Propyl Methyl Cellulose)

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1. HPMC (Hydroxy Propyl Methyl Cellulose)

Examples of film coating with HPMC

1) Example of coating using organic solution

Laboratory scale apparatus

a) Composition of coating solution

PHARMACOAT 606	8 parts
Ethanol	75.2 parts
Water	18.8 parts

b) Coating conditions

Apparatus	New Hi-Coater HCT-48N (Freund corporation)
Dimension	480 mm
Pan speed	16 (20 at later atage) min ⁻¹
Spray gun	air spray gun × 1 (AT type nozzle diameter 1.2 mm)
Drying air	55 ~ 60°C
Air flow rate	3 m ³ /min.
Spray speed	40 g/min.
Nozzle air	120 L/min.
Pressure of spray air	350 kPa
Tablet bed air temperature	46 - 39 °C
Charge per batch	25,000 tablets, 5 kg (dosage form; 8 mm, 200 mg/tab.)

c) Results

Coating time	62.5 min.
Coating solution consumption	2500 g
PHARMACOAT 606 consumption	6mg/Tab.
Disintegration time (USP disintegration test, average)	
	before cooaing; 2 min 20 sec.
	after coating; 3 min 30 sec.

Production scale apparatus

a) Composition of coating solution

PHARMACOAT 606	6 parts
Ethanol	75.2 parts
Water	18.8 parts

b) Coating conditions

Apparatus	New Hi-Coater HCT-130N (Freund corporation)
Dimension	1300 mm
Pan speed	5 (8 at later atage) min ⁻¹
Spray gun	air spray gun × 3 (AT type nozzle diameter 1.2 mm)
Drying air	50 ~ 60°C
Air flow rate	10 ~ 15 m ³ /min.
Spray speed	140 g/min.×3
Nozzle air	140 ~ 150 L/min.
Pressure of spray air	400 kPa
Tablet bed air temperature	37 - 30 °C
Charge per batch	600,000 Tablets, 120 kg (dosage form; 8 mm, 200 mg/tab.)

c) Results

Coating time	143 min.
Coating solution consumption	60.0 kg
PHARMACOAT 606 consumption	6 mg/Tab.
Disintegration time (USP disintegration test, average)	
	before cooaing; 2 min 20 sec.
	after coating; 3 min 30 sec.

2) Example of coating using aqueous solution

Laboratory scale apparatus

a) Composition of coating solution

	606	645
PHARMACOAT	6 parts	10 parts
Water	94 parts	90 parts

b) Coating conditions

Apparatus	New Hi-Coater HCT-48N (Freund corporation)	
Dimension	480 mm	
Pan speed	16 min ⁻¹	
Spray gun	air spray gun X 1 (AFT type nozzle diameter 1.2 mm)	
Drying air	70°C	←
Air flow rate	2.5 m ³ /min	←
Spray speed	20 g/min	←
Nozzle air	150 L/min	←
Pressure of spray air	200 kPa	←
Tablet bed temperature	39°C	←
Charge per batch	5 kg	←
Dosage form	6.5 mm, 120 mg/T	←

c) Results

Coating time (3% coating based on tablet weight)	83 min	50 min
Coating solution consumption	2490 g	1500 g
PHARMACOAT 606 consumption	3.6 mg/tab.	←
Disintegration time (USP disintegration test, average)		
	before coating; 2 min 20 sec	←
	after coating; 3 min 30 sec	←

Production scale apparatus

a) Composition of coating solution

	606	645
PHARMACOAT	6 parts	10 parts
Water	94 parts	90 parts

b) Coating conditions

Apparatus	New Hi-Coater HC-130N (Freund corporation)	
Dimension	1300 mm	
Pan speed	8 min ⁻¹	
Spray gun	air spray gun X 3 (AT type nozzle diameter 1.2 mm)	
Drying air	80°C	←
Air flow rate	15 m ³ /min	←
Spray speed	80 g/min×3	←
Nozzle air	170 L/min	←
Nozzle air + pattern air	250 L/min	←
Tablet bed temperature	46°C	←
Charge per batch	120 kg	←
Dosage form	6.5 mm, 120 mg/T	←

c) Results

Coating time (3% coating based on tablet weight)	286 min	171 min
Coating solution consumption	60.0 kg	35.9 kg
PHARMACOAT 606 consumption	3.6 mg/tab.	←
Disintegration time (USP disintegration test, average)		
	before coating; 2 min 35 sec	←
	after coating; 3 min 40 sec	←

1. HPMC (Hydroxy Propyl Methyl Cellulose)

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
Fig. 12: The results of dissolution tests of coated tablets at various pH values



Fig. 11: The 70% dissolution time of tablets coated with various coating materials

1. HPMC (Hydroxy Propyl Methyl Cellulose)

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*Dissolved in granulating liquid

Conditions (fluidized bed)

		Dry mixing	Solution
Granulation	Machine	Flow coater FLO-5 (Freund Corp.)	
	Drying air temperature	82-83 °C	83-84 °C
	Exhaust air temperature	28-30 °C	29-32 °C
	Air flow	2.2 m ³ /min	40 m ³ /min
	Spray gun	Schlick nozzle 1.8 mm	
	Spray air pressure	300 kPa	
	Gun position	25 cm	40 cm
	Spraying time	25 min	20 min
	Shaking/interval	6 sec/30 sec	
	End point	When the outlet temperature came to 35°C	
Drying	Drying air temperature	The same as at granulation condition	
	Air flow	2.2 m ³ /min	
	Shaking/interval	6 sec/180 sec	
	Drying time	15 min	7 min

Cumulative particle size distribution